

INTRAGASTRIC BALLOON

Corporea



Lose weight
without feeling
hungry!



A quick alternative for weight loss

AN EASY AND FASTER Method for weight loss

**It will bring direct
benefits to your
health and well-being,**
such as reduced cholesterol
and diabetes levels,
enabling the adoption of a
new food routine.



How does it work?



The Intragastric Balloon was developed to take up room in the stomach, providing a feeling of early satiety, assisting in the control of the volume of food eaten, leading to decreased caloric content, and enabling the loss of weight.



The procedure is conducted by using a digestive endoscopy.




No hospitalization is required, enabling the patient to immediately return to their daily activities, as well as there are no restrictions to the practice of physical exercises.



The removal procedure is also conducted through endoscopy at any time if the patient is not adjusted to it, without requiring any surgeries.



A multidisciplinary team is required to follow up during the treatment.



Indicated to:



Treats obesity in patients with restrictions to surgical procedures, and may result in greater weight loss in 12 months beyond expected for the 6-month treatment.¹



Previous treatment in patients with higher BMI or those who are deemed to be at high surgical risk when undergoing bariatric surgery to lose weight.



Approved for treatment of up to 12 months.

Product Characteristics:

- Manufactured with silicone, the Intragastric Balloon Corporea 12 months can be placed and removed through endoscopy, is filled with 400mL to 700mL of saline solution and methylene blue, and can remain for up to 12 months.
- Manufactured with biocompatible raw material to be used in long-term implants.
- Manufactured following a strict quality control, meeting the requirements of Good Manufacturing Practices, which assures an excellent performance to the product.





Access our website and social media

 www.balaocorporea.com.br

 [/balaocorporea](https://www.facebook.com/balaocorporea)

 [@balaocorporea](https://www.instagram.com/balaocorporea)

MEDICONE INTRAGASTRIC BALLOON - ANVISA REGISTRATION #: 80020550068
IMPORTANT: AVAILABLE BY MEDICAL PRESCRIPTION ONLY

CONTRAINDICATIONS:

Use of the Intragastric Balloon is contraindicated in the following cases:

- Pregnant women or women during the breastfeeding period;
- Handicapped individuals or unwilling to meet the restrictions associated with the diet or with the medical monitoring and their respective instructions during the program to use the Intragastric Balloon;
- Patients suffering from psychical disturbances;
- Patients younger than 18 years of age;
- Patients suffering from severe kidney and/or liver disorders;
- Patients submitted to previous gastric or previous intestinal surgery;
- Patients suffering from inflammatory disorders of the gastrointestinal tract, gastric ulcer, duodenal ulcer or specific inflammations, such as Crohn's disease, or with disposition to gastrointestinal bleeding in the upper tract, such as esophageal or gastric varicose veins, or acquired intestinal telangiectasia;

- Individuals suffering from severe cardiopulmonary or organic diseases;
- Individuals suffering from congenital or acquired anomalies of the gastrointestinal tract, such as atresia or stenoses; with large hiatal hernia;
- General chemically dependents;
- Patients suffering from any type of infection in their body;
- HIV-infected patients;
- Patients with drug administration that could cause any type of irritation or gastric complication;
- Patients with reflux esophagitis;
- Patients with stenosis or esophageal diverticulum;
- Patients with potentially hemorrhagic injuries, such as varicose veins and;
- Patients using any anti-inflammatory drugs, anticoagulants, alcoholic beverages or drugs

SINGLE USE AND NON-REUSABLE.

Reference: 1. Brooks, J. ; Srivastava, E. D.; Mathus-Vliegen, E. M. H. One-year Adjustable Intragastric Balloons: Results in 73 Consecutive Patients in the UK. Springer Science+Business Media New York: Springer, 2014. 7 p. Disponível em: <https://link.springer.com/article/10.1007/s11695-014-1176-3>.

www.medicone.com.br

Manufactured by: Medicone Projetos e Soluções para a Indústria e a Saúde Ltda.
Av. das Indústrias, nº 1585 | Distrito Industrial | Cachoeirinha - RS | CEP: 94930-230
sac@medicone.com.br | 0800 722 2728 | Phone/Fax: +55 (51) 3470-0800

Code/Review: 67381/00

